



# Aspect Analytics

Dedicated Software to Support High Throughput  
Spatial Multi-Omics Applications

Marc Claesen, PhD  
Co-founder and Chief Executive Officer  
[marc.claesen@aspect-analytics.com](mailto:marc.claesen@aspect-analytics.com)

AGENTSCHAP  
INNOVEREN &  
ONDERNEMEN



Vlaanderen  
is ondernemen



Accelerated by  
Health

 **NVIDIA.**  
INCEPTION PROGRAM



  
THE MICHAEL J. FOX FOUNDATION  
FOR PARKINSON'S RESEARCH



# Presentation outline

1. Introduction
2. Our spatial multi-omics approach
3. Case studies
4. Conclusion

# Aspect Analytics in a nutshell

## Bioinformatics solutions provider for spatial biology

- **Weave cloud platform** for spatial biology data management and analysis
  - **Fit-for-purpose tools**: customizable to fit core lab workflows
- **Bioinformatics services**: multi-omics integrations, data-driven QC, software co-development, ...



## Our mission


- Provide a **unified platform solution** for spatial biology data analysis
- Unlock full potential of **spatial multi-omics in high-throughput contexts**



# Spatial multi-omics

Improving our understanding of tissue biology

- **Spatial context is a key driver** of complex biology
- **Multi-omics needed** to capture full biological process






# Spatial multi-omics

Improving our understanding of tissue biology

- **Spatial context is a key driver** of complex biology
- **Multi-omics needed** to capture full biological process

Data analysis challenges in spatial multi-omics

- **Multiple assays needed to target all omics layers**
  - Integrate data across sections and vendor formats
- **Large data volumes and complexity**
  - Multi-petabyte storage & on-demand computations
- **Interdisciplinary collaboration** with many stakeholders
  - Team communication drives data analysis efficiency



# Broad assay coverage

## Spatial transcriptomics

- **Subcellular assays:** Xenium, CosMx, MERSCOPE
- **Minibulk assays:** Visium, Visium HD, GeoMx, Seeker

## Spatial proteomics

- **mIF based:** COMET, PhenoCycler, CellDIVE, CellScape
- **Other approaches:** Hyperion (IMC), MALDI-IHC

## Mass spectrometry imaging

- **instruments:** timsTOF fleX, solariX, scimaX, MRT, Orbitrap, Synapt, iMScope (+ imzML)
- **analytes:** drugs, lipids, glycans, metabolites, peptides, proteins



# Broad assay coverage

## Spatial transcriptomics

- **Subcellular assays:** Xenium, CosMx, MERSCOPE
- **Minibulk assays:** Visium, Visium HD, GeoMx, Seeker

## Spatial proteomics

- **mIF based:** COMET, PhenoCycler, CellDIVE, CellScape
- **Other approaches:** Hyperion (IMC), MALDI-IHC

## Mass spectrometry imaging

- **instruments:** timsTOF fleX, solariX, scimaX, MRT, Orbitrap, Synapt, iMScope (+ imzML)
- **analytes:** drugs, lipids, glycans, metabolites, peptides, proteins

Spatial multi-omics  
data integration  
across any combo  
of these assays



Check booth #3  
for live demos!





# Our spatial multi-omics approach



# Spatial multi-omics data comes in *stacks*



# Spatial multi-omics data comes in *stacks*





Data integration is challenging due to differences in

- **Coordinate systems** due to spatial resolution, different sections, probe size, ...
- **Data formats**: each assay has its own proprietary vendor format





# weave platform for spatial biology





morphology	pathologist annotations	spatial transcriptomics	spatial proteomics	spatial metabolomics	other data
_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____

## Stack fusion

- Integrate readouts across assays and sections to obtain unified, holistic dataset
- Enables true spatial multi-omics using any number and combination of assays



# weave platform for spatial biology




morphology	pathologist annotations	spatial transcriptomics	spatial proteomics	spatial metabolomics	other data
_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____

## Stack fusion

- Integrate readouts across assays and sections to obtain unified, holistic dataset
- Enables true spatial multi-omics using any number and combination of assays
- Facilitate downstream analysis and interactive visualization via dataframe-like interfaces



multimodal visualization





unified data structure













- Cell Seg...  
 Opacity:
- MALDI ...  
 Opacity:
- H&E (Ma...  
 Opacity:
- H&E Mo...  
 Opacity:
- H&F (M...










# Cloud platform for spatial biology analysis



# Case studies

# Spatial profiling of prostate cancer (PCa)

Goal: improving understanding of PCa biology via spatial biology

- Roughly **1 in 8 men will be diagnosed with PCa** during their lifetime
- **PCa is a very complex disease**
  - high inter -and intra-tumoral heterogeneity
  - multiple morphological disease states in single section
  - involves molecular changes across multiple omics classes



Zhang et al., bioRxiv 2023.08.28.55505  
<https://doi.org/10.1101/2023.08.28.555056>

Wanqiu Zhang  
Xander Spotbeen  
Sebastian Vanuytven  
Sam Kint  
Tassiani Sarretto  
Fabio Socciarelli  
Katy Vandereyken  
Jonas Dehairs  
Jakub Idkowiak  
David Wouters  
Jose Ignacio A. Larizgoitia  
Gabriele Partel  
Alice Ly  
Vincent de Laat  
Maria José Q Mantas  
Thomas Gevaert

Wout Devlies  
Chui Yan Mah  
Lisa M Butler  
Massimo Loda  
Steven Joniau  
Bart De Moor  
Alejandro Sifrim  
Shane R. Ellis  
Thierry Voet  
Marc Claesen  
Nico Verbeeck  
Johannes V. Swinnen



# Spatial profiling of prostate cancer (PCa)

Goal: improving understanding of PCa biology via spatial biology

- Roughly **1 in 8 men will be diagnosed with PCa** during their lifetime
- **PCa is a very complex disease**
  - high inter -and intra-tumoral heterogeneity
  - multiple morphological disease states in single section
  - involves molecular changes across multiple omics classes



KU Leuven  
Institute for  
Single Cell  
Omics



UOW  
AUSTRALIA




THE UNIVERSITY  
of ADELAIDE

Zhang et al., bioRxiv 2023.08.28.55505  
<https://doi.org/10.1101/2023.08.28.555056>

## High-level design:

- 16 samples coming from 8 patients - 8 malignant and 8 non-malignant biopsies
- Histology and multi-omics combination:
  - Spatial analysis: **spatial transcriptomics** (Visium) + **spatial lipidomics** (Orbitrap)
  - Bulk: **snRNAseq** + **lipidomics** (LC-MS)


# Stack design and experimental approach



# Step 1: multimodal registration pipeline

## Approach


- **Inter-section registration via H&E's**
  - H&E available for both Visium & MSI sections
  - Shared modality facilitates registration
- **Registrations per section of omics and H&E**
  - MSI to MSI H&E, Visium to Visium H&E
- All **registrations are non-rigid** to account for serial sections and complex sample deformation



## Outcome

- **Full coordinate transformation available from between any layer combination in the stack**





**Visium Assay**  
Opacity: [Slider]

**H&E (MALDI MSI Tissue Section)**  
Opacity: [Slider]

**H&E (Visium Tissue Section)**  
Opacity: [Slider]

**MALDI Mass Spectrometry Imaging of Lipids**  
Colormap: magma  
Domain: Min/Max  
Opacity: [Slider]  
Zero Transparent:

PE 40:3 [M+Na]<sup>+</sup> - m/z 820.5836 - ppm  $\epsilon$   [Slider]

PE O-34:2 / PE P-34:1 [M+H]<sup>+</sup> - m/z 702.  [Slider]

PE 42:7 [M+H]<sup>+</sup> - m/z 818.5716 - ppm  $\epsilon$   [Slider]

Search

- AGPAT5
- AGPS
- AGR2
- AGR3
- AGRN
- AGT
- AGTPBP1
- AGTR1
- AGTRAP
- AHCTF1
- AHCY
- AHCYL1

Cell Sets

- Visium Spot Labelling




In this project a collection of human prostate samples has been analyzed using spatial multi-omics, namely spatial lipidomics via MALDI MSI (Thermo), spatial transcriptomics using Visium (10x Genomics) and traditional histology.

## Step 2: data integration to Visium spots


### Principled strategy for data integration

- **Different spatial footprints** between assays
  - MSI: 30 micrometer per pixel
  - Visium: 55 micrometer per spot
- We matched MSI pixels to Visium spots by **weighing their spatial overlap**



### Outcome

- **Integrated data structure for downstream analysis**
  - Aggregated, representative mass spectrum for each Visium spot



Transcripts	Lipids












# Spatial correlation / anticorrelation

Correlation: 0.694

PC 34:2 [M+H]<sup>+</sup>



TMEFF2




Correlation: -0.631

SM 34:1;O2 [M+H]<sup>+</sup>


AGR2








# Lipid: PE 40:2 [M+H]<sup>+</sup>




# Gene: MSMB




High positive correlation in normal samples

NORMAL SAMPLES

# Lipid: PE 40:2 [M+H]<sup>+</sup>



# Gene: MSMB



High positive correlation in normal samples


NORMAL SAMPLES




Lipid: PE 40:2 [M+H]<sup>+</sup>

Gene: MSMB

NORMAL  
SAMPLES




High positive  
correlation in normal  
samples



CANCER  
SAMPLES



Almost no correlation  
in cancer samples






Lipid: PE 40:2 [M+H]<sup>+</sup>

Gene: MSMB



NORMAL  
SAMPLES

CANCER  
SAMPLES



High positive  
correlation in normal  
samples


Almost no correlation  
in cancer samples




Lipid: PE 40:2 [M+H]<sup>+</sup>

Gene: MSMB


NORMAL  
SAMPLES




High positive  
correlation in normal  
samples



CANCER  
SAMPLES



Almost no correlation  
in cancer samples





# High-grade serous ovarian cancer (HGSC)

Ovarian cancer is the most lethal gynecologic malignancy


- HGSC is the most common histotype, accounting for over 70% of all ovarian cancer deaths

Current treatment regime:

- various combinations of surgery and chemotherapy

MD Anderson uses follow up laparoscopy (6-8 weeks) to look for minimum residual disease (MRD)

- **No biomarkers to predict development of MRD or mechanism of chemoresistance**



# Materials + Methods - Tissue stack generation

## Cohort (8 samples)

- Eight samples
  - Four MRD+
  - Four MRD-
- FFPE
  - 10 $\mu$ m serial sections
  - All sections for analysis collected at same time ("stack")

## Measurement Stack (6 modalities)



10x Visium

### Visium (~55 $\mu$ m)

- + H&E
- Spatial Transcriptomics

THE UNIVERSITY OF TEXAS  
MD Anderson  
~~Cancer Center~~  
Making Cancer History®



Lunaphore COMET

### Multiplex IF (<1 $\mu$ m)

- 20 markers
- + DAPI

THE UNIVERSITY OF TEXAS  
MD Anderson  
~~Cancer Center~~  
Making Cancer History®



Bruker timsTOF fleX

### MSI (20 $\mu$ m) - Erin Seeley

- metabolites
- glycans
- peptides
- post -measurement H&E

 **TEXAS**  
The University of Texas at Austin




# Why combine modalities?


**Immunofluorescence (IF)** provides high resolution images where we can **identify** and spatially segment **individual cells**

The IF markers allow us to find the **“type”** of the cell using **previous knowledge**



**MSI** provides **novel molecular information** that can only be obtained using the technique



Describe novel spatial environments by molecular signals





# Step 1: multimodal co-registrations







# Step 1: multimodal co-registrations








# Step 1: multimodal co-registrations








# Step 2: matching MSI pixel footprint



MSI  
pixels  
1-n...

*Aggregation to the MSI pixel level...*

cell type proportion	antibody markers	MSI peptides	MSI metabolites	MSI glycans
				



weave











# Differential MSI signals by cell type aggregation (tumor vs CAFs)

■ Tumor

m/z 1742.5924 -  
Hex5HexNAc4NeuGc1\_mouse

m/z 1325.757 - histone H4

m/z 128.0354 - Pyroglutamic acid




■ Tumor


■ Cancer Associated Fibroblasts

■ Immune cells

logfc ~0.53




logfc ~0.92




logfc ~0.69



# Spatial exploration of Tumor-Stroma interface




-1000 µm



Interface

+1000 µm





# Xenium reproducibility analysis


Goal: assess reproducibility of Xenium assay

- Analyzed **Xenium data acquired from 2 serial sections**
  - Experiments run by different operators
- Performed the analysis on **3 human carcinoma FFPE tissue types**

Breast tissue

Colon tissue

Lung tissue







# Xenium reproducibility analysis

## Leveraging Weave


- **Automatic registration between each pair of Xenium data sets**
  - Per section, DAPI to H&E stained microscopy
  - Between sections, H&E to H&E
- **Joint data visualization of each pair into a single view**



ABCC11 for each section (red/green)



MYLK transcripts per section



MYLK density per section



# Statistical check


Analysis 2: assess whether transcript count distributions are comparable

- Stratified across gross morphological regions defined by expert pathologists
- Non-parametric KS test for distribution similarity per region x gene pair


Tissues	Tumor region	Stroma region	Immune region
Breast	90%	90%	89%
Colon	98%	94%	99%
Lung	94%	100%	99%

Percentage of genes with equivalent distribution per region.





Breast tissue



Colon tissue




Lung tissue




# Xenium - Visium comparison

Analysis 3: check concordance of Xenium & Visium on same section


- Investigate gene expression in same area (here CXCL9)




CXCL9 Xenium (spots)



CXCL9 Xenium (density)



CXCL9 Visium (spots)



CXCL9 overlay



# Conclusion

# Conclusion

## Cloud platform for spatial multi-omics

- Foundation for data management and analysis
- Customized fit-for-purpose data workflows
- Vendor neutral

## Collaboration & services

- Research partner for high-end bioinformatics
- Development of new assays and applications
- Deep integration with vendors

## Addressing spatial multi-omics challenges

- Integration across analytical platforms
- Supporting high-throughput workflows
- Effective communication of spatial omics insights







# Aspect Analytics

info@aspect-analytics.com  
<https://www.aspect-analytics.com>